



1200W Diode Laser

Faster, stronger, more comfortable


For all skin types and hair colors

Permanent hair removal



Updated Treatment Handle


- Equipped with 3 removable treatment heads to meet the needs of different parts of the body
- Smart touch screen, easy to operate



Larger spot size, Higher laser power

FOCUSLIGHT laser device, high durability

Long service life



Spot size
12mm * 24mm




Unique Features

- The gold standard for permanent hair removal
- High output power **1200W**
- Equipped with special skin cooling plate (optional)
- Upgraded treatment handle with **3** removable treatment heads and smart touch screen
- Double filters to keep water quality
- 808nm operating system or 755+808+940+1064nm system optional
- **FOCUSLIGHT** best-quality laser bars
- Fast hair removal of all skin types (I, II, III, IV, V, VI) and hair colors
- Dual Cooling systems: TEC/Semiconductor + Sapphire contact cooling, bring a painless and comfortable treatment experience



Clinical Results:

Clinical tests show that in 95% of cases, a hair reduction of more than 80% is reported after 3 ~ 5 sessions. After 1 ~ 2 sessions, absolutely cases show a reduction of between 30% and 50%.



Technical Parameters:

Output power	1200W
Laser type	Diode laser
Spot size	Φ5mm, 12*12mm, 12*24mm (optional)
Frequency	1Hz ~ 10Hz
Fluence	0 ~ 120J/cm ²
Pulse duration	10ms ~ 200ms
Laser intensity	5% ~ 100%
Operating system	808nm system or 755+808+940+1064nm system optional
Cooling system	Water+Wind+sapphire+TEC/Semiconductor
Voltage	220V/110V, 50/60HZ

