


# Hydra Queen Non-invasive mesotherapy

Brilliant and bright, deep into the pores to take away the dirt, sweep the net

**NEW!**



- Nanocrystals use nanoscale contacts that can be as thin as 250 nanometers and as short as 0.1 millimeters.
- Just open the stratum corneum of the skin without touching the dermis layer.
- It not only solves the problem of skin absorption, but also makes it painless and non-invasive!





### Natural application and absorption

Natural application and absorption only penetrate into the epidermal layer. It is difficult for skin care product molecules to penetrate and the absorption amount is very low.

### Nanocrystalline nanochips

(the product can evenly penetrate into the dermal layer)

The use of nanochips only penetrates the epidermal barrier layer without touching blood vessels and nerves. Painless and non-invasive absorption is greatly improved.



### EMS function

Microelectronic thermal conductivity technology opens pores and increases product absorption rate by more than 40%



sagging skin



Dry and water shortage



Fine water wrinkles



gloomy

# Hydra Queen eight in one superbubble

Dry and wet separation  
(dual cleaning system)

Bubble Cleaning

Skin Scrubber vibration cleaning


EMS Nano mesotherapy

Hot and cold +RF

Bipolar HR

Ultrasonic lead-in

Microdermabrasion



# Hydra Queen Eight in one Device



## **Bubble Cleaning**

A new generation of instruments, breaking through the technical level of traditional small bubble instruments, can open pores, soften keratin, make skin cleaner and deeper improve great of the dark skin and blackheads, and make skin clean, white and tender.



## **EMS Nano mesotherapy**

Nanocrystalline is promoted by opening absorption channels in the skin. Absorption into skin care products, so use nanocrystalline import. When the skin care products can be more absorbed into the skin.



## **Bipolar HR**

Use for body heat therapy to make skin tissue suck received a lot of heat combined with new collagen, smooth tighten the skin.



## **Hot and cold +RF**

Cold and hot lifting, can tighten and removal wrinkle, strengthen cell regeneration, let cells recovery activity, repair sensitive skin.



## **Skin Scrubber vibration cleaning**

Whitening, remove dead skin, blackhead, Import function.



## **Ultrasonic lead-in**

Ultrasonic vibrations of 1 to 3 million per minute frequency, open the pores, make the ingredients fast rapid penetration into the basal layer of the skin, allowing absorption rate more than 90%.



## **Microdermabrasion**



## High pressure water oxygen pen

450m/s speed, 50-85kpa injection pressure, directly inject nutrients into the dermis, stimulate cell regeneration, deep hydration, lasting effect. Unique shape design, needle essence import, easy to exchange freely according to the situation, reduce the cost of consumables.


### **OXYGEN JELTION**

High pressure oxygen injection, particle penetradeep into the muscle The skin layer enhances waterretention by narrowing pores more Easy to absorb



# What's special about it?

Triple-efficiency operating head



Three major areas with clear division of labor can achieve three times the effect at one time

Vacuum negative pressure, combined with negative pressure suction, takes away all dirt from the face



# System Specification

## Functions

## Specifications

**Voltage:** 100v-120v 60Hz 200v-240v 50Hz

**Power:** 200W

**HR frequency:** 0.2-0.5 Mhz

**coldand hotscope:** Cold 0-6C Hot36-46° C

**Ultrasound:** 1.1MHz

**Probe:** 3 soft heads,8 hard heads

**Clean water mill:** Positive pressure 4-5KgNegative pressure 96Kpa.Flows50L

**Display:** 10.1inches Resolution:1024X600

**Solution bottle:** A.Deep layer cleaning  
B.Whitening and moisturizing  
C.Deep cleaning  
D.Instrument piping cleaning

**Skin scrubber vibration frequency:** 26-30MHz

**Ems Beauty plastic head size:** 4X4mmX250Um

**Roller Speed up:** frequency200Hz microcurrent50-800(Ua)Display screen:10.1inches

**Visual probe:** dry head 3 (large, medium, small)

# Before —After



Before    .....    After

Before    .....    After

.....



Before    .....    After

Before    .....    After

Before    .....    After



Before    .....    After

Before    .....    After

Before    .....    After